



# High-Quality Collagen Solutions

## Solution type

High-Quality Collagen excellent for all types of research including cell biology, biochemistry, biophysics, tissue engineering, etc.

0.45 µm filtered



### Acid Soluble Collagen

- Strong gel strength
- Suitable for 3D-cell culture
- Superior fibril formation with remaining telopeptides and crosslinks.

### Pepsin Solubilized Collagen

- Both N- & C-telopeptides are removed
- Suitable for coating culture vessels

### Type I Collagen

- Most abundant collagen in vertebrates
- Major component in collagen fibrils
- More than 95% pure

### Type III Collagen


- Abundant in skin tissues
- Rich in juvenile tissues and is important in wound healing
- Non-collagenous domains are removed by pepsin treatment
- P#892 107&108 contain about 20% of type I collagen

### Type V Collagen


- Fibrillar collagen
- Rich in the cornea
- Relatively rich near the basement membrane
- Abundant in fine collagen fibrils
- Non-collagenous domains are removed by pepsin treatment

#### Usage Examples


##### Collagen coating for a 24 well plate



##### Collagen gel prep for a 24 well plate





## Reconstituted fibrils



Type I collagen fibrils reconstituted at 37°C under physiological conditions. SEM image (5000x)

## Excellent cell adhesive property



Adhesion of human keratinocyte cell line FEPE1L-8 at 1 day



# Collagen for research reagents

## Powder type

### Features of Collagen Powder

- Easy to adjust the concentration
- Able to use various solvents
- Native triple helical structure retained

### How to dissolve

- Dissolve in an acidic solution, such as 5 mM acetic acid or 1 mM hydrochloric acid.
- To increase the concentration of collagen, first dissolve in water, and then add the required amount of acid solution.
- When the concentration exceeds 5 mg/mL, handling becomes difficult due to the high viscosity.
- Can be dissolved up to a concentration of 10 mg/mL

### Applications

- 2D and 3D cell culture substrates
- Bio-ink
- Production of collagen moldings
- R&D for drug delivery system (DDS)



Product code	Product name	Capacity
<b>Solution type</b>		
892 101	Type I collagen, Bovine skin, Acid soluble, 3mg/mL	20mL
892 102		100mL
892 103	Type I collagen, Bovine skin, Pepsin-solubilized, 3mg/mL	20mL
892 104		100mL
892 107	Type III collagen, Bovine skin, Pepsin-solubilized, 3mg/mL	5mL
892 108		20mL
892 151	Type V collagen, Bovine cornea, Pepsin-solubilized, 3mg/mL	1mL
892 111	Type I collagen, Porcine skin, Pepsin-solubilized, 3mg/mL	20mL
892 112		100mL
<b>Powder type</b>		
892 140	Type I collagen, Bovine skin, Acid soluble	100mg
892 141		500mg
892 142	Type I collagen, Bovine skin, Pepsin-solubilized	100mg
892 143		500mg
892 144	Type I collagen, Porcine skin, Pepsin-solubilized	100mg
892 145		500mg

Manufactured by

**nippi** Nippi, Incorporated

1-1-1 Senjumidori-cho, Adachi-ku, Tokyo 120-8601, JAPAN

PHONE : +81-3-3888-5184

FACSIMILE : +81-3-3888-5136

E-mail : protein-info@nippi-inc.co.jp

Agency: